



SleepWeaver®

Soft Cloth Masks



Anew®



Élan®



3D



Advance



Advance
Small



Advance Pediatric

Soft, Quiet, Comfort!

SleepWeaver is the world's first line of all-cloth Positive Airway Pressure (PAP) masks. Created for CPAP users by CPAP users.

- **CLOTH** Inflates like a balloon and conforms to your facial features.
- **CLOTH** is soft, light, quiet, and does not cause facial irritation from plastic.
- **CLOTH** creates a comfortable loose seal that promotes your compliance to CPAP Therapy.
- **CLOTH** is as comfortable as your favorite pajamas.
- **CLOTH** moves with you, so you can sleep more naturally.
- **CLOTH** conforms and seals perfectly to your face for more effective air delivery.
- **CLOTH** wicks away moisture to let your skin breathe.
- No buckles, clips, or faceplates touching and irritating the skin.
- No more noisy air blowing into your bed partner's face.
- Easy to clean and maintain!

Approved for reimbursement by major medical insurance plans including Medicare.

*Soft **CLOTH** is Better!*

SleepWeaver®

Advance

Soft Cloth Masks

Our flagship product



- Zzzephyr Seal® conforms to your face to provide a soft and comfortable seal.
- 3-Point soft headgear.
- Design allows for wide range of adjustability.
- One size fits most.
- Accommodates eyeglasses – great for reading or watching TV!



Choose a color and make it your own!

SleepWeaver®

Advance Small

Soft Cloth Masks

For smaller facial features



- 15% smaller than our SleepWeaver Advance Mask.
- Fits petite faces.
- Enhanced seal offers improved comfort.
- 3-Point soft headgear.
- Design allows for wide range of adjustability.



Available In Beige.

SleepWeaver®

Advance Pediatric | Soft Cloth Masks

The world's first and only soft cloth pediatric PAP mask.



- FDA cleared for patients between 2-7 years old.
- 5-point adjustable headgear.
- Constructed of soft cloth that wicks away moisture.
- Ideal for sensitive skin.
- Enhanced seal for a softer feel on your face.



Available In Blue.

SleepWeaver®

Élan® | Soft Cloth Masks

The most adjustable mask in our lineup



- Quick, simple, easy, and cost effective mask cushion replacement.
- 90° swivel connector with 360° rotation allows for more natural sleeping positions.
- Open-face head strap design accommodates eyeglasses.
- Tether strap compatible to position tube overhead.



Choose a color and show off your style!

Available in Large, Regular & Small

SleepWeaver®

3D | Soft Cloth Masks

The newest soft option



Available In Blue.

Available in
Large & Regular

- The flexible elastic interface material allows for sealing on many face types.
- Wing-shaped stay routes the headgear under the ears – ensuring comfort and a great fit.
- Rounded edge straps – reduces strap marks.
- Permanent 360° swivel provides reliable performance for the tubing in any position.

SleepWeaver®

Anew® | Soft Cloth Masks

The only full-face soft cloth PAP mask in the world!



Available In Blue.

Available in Large,
Regular & Small

- Dual interface design securely seals around both the nose and the mouth. Perfect for mouth breathers.
- 90° swivel connector with 360° rotation allows for more natural sleeping positions.
- Exhalation holes in the front ensure a quiet and gentle flow of air.
- Super lightweight and easy-to-clean design.

SleepWeaver®

Accessories



18"
Feather
Weight
Tube®

Hose Clip



Note: SleepWeaver Feather Weight Tubes are designed to fit any PAP mask with a 22 millimeter connection.



101653 REV 3



Contact us for more information about the SleepWeaver line of Soft Cloth PAP Masks:

1-888-825-9640 • 1-724-858-2837
info@circadiance.com • www.circadiance.com