

Deep restorative sleep for mouth breathers.

The World's First Soft Cloth Full Face PAP Mask



Comfortable Seal for Mouth Breathers:

The dual interface design individually seals around the nose and the mouth allowing only soft cloth to touch the skin.

Open Face Headgear Design:

Allows virtually unrestricted vision; great for people who like to fall asleep while reading or watching TV. Accommodates eyeglasses.

90° Elbow Swivel Connector with 360° Rotation:

Allows users to sleep in more natural positions; easy to remove and clean.

Gentle Exhalation:

Exhalation holes are in the front of the mask providing a quiet and gentle flow of air.

Easy to Clean:

Quick & simple to clean the material that routinely comes in contact with the skin surface.

Super Lightweight:

Users barely notice the weight of this mask.

3 Sizes:


Multiple size options allow users to find the perfect fit for their face; expands the number of users the mask will fit.

Moisture Vapor Breathable:












High-tech material allows the skin to breathe naturally, wicks away moisture and leaves no pressure marks that are common with rubber or silicone seals.

Loose Fit:

As with all SleepWeaver masks, the unique "loose fit" provides ultimate comfort and eliminates strap marks.



SleepWeaver Anew Ordering Information

Product	Photos	Sizes	HCPCS Codes	Part Numbers	Color Options		
SleepWeaver Anew Mask & Headgear		Small Regular Large	A7027 - Mask A7035 - Headgear	Small	100955	Blue	
				Regular	100960	Blue	
				Large	100965	Blue	
SleepWeaver Anew Mask Only		Small Regular Large	A7027 - Mask	Small	100970	Blue	
				Regular	100975	Blue	
				Large	100980	Blue	
SleepWeaver Anew Headgear Only		Small Regular	A7035 - Headgear	Small	101145	Black	
				Regular	101084	Black	

SleepWeaver Anew Accessories

					
Tether Strap - 100638	Hose Clip - 100775	DVD 10 Pack - 101086	18 inch Tube - 100477	Washer - 100778	Connector - 100777
Tether Strap attaches to Tube & Headgear	Hose Clip prevents Tube movement	SleepWeaver Anew Instructional DVD video with printable Instructions for Use	SleepWeaver Feather Weight Tube	SleepWeaver Anew Washer Replacement Part	SleepWeaver Anew Threaded Connector Replacement Part

Contact us for more information about the SleepWeaver line of Soft Cloth PAP Masks:

1-888-825-9640 * 1-724-858-2837

orders@circadiance.com

www.Circadiance.com



SleepWeaver
Soft Cloth PAP Masks



101133 REV3