

The Best in Comfort and Technology



- Soft
- Light
- Quiet
- Easy to fit
- No skin irritation from plastic
- Comfortable seal on many face types
- Compatible with common CPAP/Bi-Level machines
- Durable
- Easy to clean
- Economical

"Love the mask"

Raymond B.

***"Excellent Mask,
Very Comfortable"***

Rich H.

***"I fell in love with this
one (SleepWeaver 3D)"***

John B.

***"Best I have used,
works very well and
comfortable."***

Mary Jane S.

A new and more flexible elastic interface material:

This creates our best seal ever on many face types.

New rounded edge straps:

Reduces strap marks.

New wing-shaped stay:

Routes the headgear under the ears, ensuring comfort and a great fit.

Permanent 360° Swivel:

Provides reliable performance for the tubing in any position, including over the headboard.

Nothing Breathes Like Cloth:

Breathable cloth design eliminates fogging common in hard plastic masks. Keeps skin dry and comfortable.

Lightweight:

What could be lighter than a thin piece of cloth?

Accommodates Eye Glasses:






For those who want to read or watch TV while in bed, SleepWeaver 3D allows users to wear eye glasses with the mask on.

The Joy of Soft Cloth:


No one can deny preferring a soft cloth against their face as opposed to rubber, silicone and other gels. Cloth masks eliminate the need for a cushion or nasal pad.



SleepWeaver 3D ordering information

Product	Photos	HCPCS Codes	Part Number
3D Mask & Large Headgear		A7034-Mask A7035-Headgear	101491
3D Mask & Regular Headgear		A7034-Mask A7035-Headgear	101486
3D Mask Only		A7034-Mask	101496
3D Large Headgear Only		A7035-Headgear	101517
3D Regular Headgear Only		A7035-Headgear	101516
Feather Weight Tube Only		A7037-Tube	100477
Hose Clip		N/A	100775

SleepWeaver 3D Sizing Guide



Watch SleepWeaver 3D Videos Online:
www.circadiance.com/3dvideos

101563 REV 1

Phone: 888-825-9640 or +1-724-858-2837
orders@circadiance.com www.circadiance.com



SleepWeaver
 Soft Cloth PAP Masks